



ZipCombo Microhematocrit Centrifuge w/ 12 Place Hematocrit Rotor

The LW Scientific ZipCombo centrifuge provides all the performance of popular “Stat” centrifuges at a fraction of the cost. This convertible centrifuge can be used for blood, urine, DNA separation, and HPLC microfiltration. Spin down hematocrits at 12,000 rpm in only 3 minutes. Optional rotor is available to spin micro-tubes at 12,000 rpms for serum in 3 minutes. Use the digital controls to decrease the speed to also spin urines in micro-tubes. Featuring a brushless motor with multiple rotor choices, along with a safety switched transparent lid. Variable speed of 0-12,000rpm with Maximum RCF=6900g. Included is a reader card and 200 Mylar-reinforced 40mm hematocrit tubes.

100-240V 50/60Hz auto-switching

One year warranty

Distributed by ServiceNet Medical X-ray 1-800-830-1035